


Tip-sample interaction on graphite studied in the thermal oscillations regime

Giovanna Malegori




Dipartimento di Matematica e Fisica, Università Cattolica del Sacro Cuore, Brescia, Italy
Dipartimento di Fisica, Università degli Studi di Milano, Milano, Italy

- Introduction
- Probing the tip-sample interaction
 - Frequency shift
 - Potential from Boltzmann distribution
 - Mean-square displacement from power spectral density

AFM cantilevers display thermal flexural oscillations

- noise source
- calibration of the cantilever spring constant $k^{[*]}$
- **probing van der Waals tip-sample interaction**



[*] J.E.Sader *et al.*, *Rev. Sci. Ins.* **70**, 3967 (1999)

- The **free AFM cantilever** oscillates due to random thermal excitations
- As the tip approaches the sample surface, the **tip-sample interaction** deflects the cantilever and modifies its thermal vibrations → tip-sample force gradient
- The **jump-to-contact** occurs due to liquid meniscus pulling the tip (humidity) [*]



[*] M. Luna *et al.* *J. Phys. Chem B* **103**, 9576 (1999)

Free cantilever thermal motion →
Langevin equation

$$m^* \frac{d^2 u}{dt^2} + \gamma \frac{du}{dt} + ku = F_{rand}$$

- u cantilever displacement
- k spring constant
- m^* effective mass
- γ damping
- F_{rand} thermal stochastic force



$$\langle F_{rand}(t) \rangle = 0$$


$$\langle F_{rand}(t_1), F_{rand}(t_2) \rangle = \delta(t_1 - t_2)$$

V.L.Mironov *Fundamentals of SPM* (2004)
 F.J.Giessibl, *Rev. Mod. Phys.* **75**, 949 (2003)
 D.T.Gillespie, *Am. J. Phys.* **61**, 1077 (1993)

Cantilever near the surface → F_{ts} tip-sample interaction force

$$m^* \frac{d^2 u}{dt^2} + \gamma \frac{du}{dt} + ku = F_{rand} + F_{ts}(z)$$

$$F_{ts}(z) = F_{ts}(z_0) + \frac{\partial F_{ts}(z_0)}{\partial z} u + \dots$$



- The constant force $F_{ts}(z_0)$ only displaces the equilibrium position z_0 (static deflection)
- The force derivative $\partial F_{ts} / \partial z$ influences the cantilever oscillations



Interacting cantilever

A.P.E. Research
NANOTECHNOLOGY

Small amplitude oscillations →
force gradient $\partial F_{ts} / \partial z$ constant during the cantilever oscillations

$$m^* \frac{d^2 u}{dt^2} + \gamma \frac{du}{dt} + ku = F_{rand} + F_{ts}(z_0) + \frac{\partial F_{ts}(z_0)}{\partial z} u$$

$$m^* \frac{d^2 u'}{dt^2} + \gamma \frac{du'}{dt} + k^* u' = F_{rand} \quad k^* = k - \frac{\partial F_{ts}}{\partial z}$$

u' cantilever displacement from the new equilibrium position z_0


The tip-sample force gradient is directly related to the effective spring constant k^*

V.L.Mironov *Fundamentals of SPM* (2004)
F.J.Giessibl, *Rev. Mod. Phys.* **75**, 949 (2003)

11 March 2010

MICROSCOPIES A SONDE LOCALE
FORUM 2010


- Introduction
- Probing the tip-sample interaction
 - Frequency shift
 - Potential from Boltzmann distribution
 - Mean-square displacement from power spectral density



[1] O.Pfeiffer *et al.*, *Phys. Rev. B* **65**, 161403 (2002)
 T.Drobek *et al.*, *Phys. Rev. B* **64**, 045401 (2001)

- [2] Y.Sugimoto *et al.*, *Appl. Phys. Lett.* **91**, 093120 (2007)
 M.A.Lantz *et al.*, *Science* **291**, 2580 (2001)
 H.Hoelscher *et al.*, *Phys. Rev. Lett.* **83**, 4780 (1999)
- [3] D.O.Koralek *et al.*, *Appl. Phys. Lett.* **76**, 2952 (2000)
 A.Roters *et al.*, *J. Phys. Condens. Matter* **8**, 7561 (1996)
- [4] J. Buchoux *et al.*, *Nanotechn.* **20**, 475701 (2009)

Experimental setup



[*] J.L.Hutter, *Langmuir* **21**, 2630 (2005)


- One temporal trace of the cantilever thermal motion
- Three different approaches to measure the tip-sample interaction
 - Frequency shift method
 - Potential from Boltzmann distribution
 - Mean-square displacement from power spectral density

- Introduction
- Free cantilever thermal vibrations
- Probing the tip-sample interaction
 - Frequency shift
 - Potential from Boltzmann distribution
 - Mean-square displacement from power spectral density



Free cantilever's flexural modes

A.P.E. Research
NANOTECHNOLOGY






Boundary conditions

$$f_n^F = \alpha_n^2 \frac{h}{L^2} \sqrt{\frac{E}{12\rho}}$$

Geometrical parameters Material properties

Mode	1	2	3	4	5
f_n (kHz)	10.8	74.3	215	431	724
f_n / f_1 (exp.)		6.9	20	40	67
f_n / f_1 (th. [*])		6.27	17.55	34.39	56.84



[*] H.J. Butt *et al.*, Nanotechnology 6, 1 (1995)

11 March 2010



Frequency shift method

The resonance frequency of free cantilever is

$$f_0 = \frac{1}{2\pi} \sqrt{\frac{k}{m^*}}$$

k spring constant, m^* effective mass

The tip-sample interaction F_{ts} changes the resonance frequency.
For small amplitude oscillations

$$k^* = k - \frac{\partial F_{ts}}{\partial z} \quad f = \frac{1}{2\pi} \sqrt{\frac{k^*}{m^*}}$$

k^* effective spring constant

The frequency shift is directly related to the tip-sample force gradient

$$f = f_0 + \Delta f$$


$$\frac{\partial F_{ts}}{\partial z} = -2k \frac{\Delta f}{f_0}$$

Y.Martin *et al.*, *J. Appl. Phys.* **61**, 4723 (1987)

F.J.Giessibl, *Rev. Mod. Phys.* **75**, 949 (2003)




HOPG







HOPG




Power
Spectral
Density






HOPG




Power
Spectral
Density





Thermal PSD approaching the HOPG surface



HOPG




Power
Spectral
Density






HOPG




Power
Spectral
Density







HOPG




Power
Spectral
Density







HOPG




Power
Spectral
Density








HOPG




Power
Spectral
Density







HOPG




Power
Spectral
Density







HOPG




Power
Spectral
Density







HOPG




Power
Spectral
Density






HOPG





The jump-to-contact occurs when

$$\left. \frac{\partial^2 V_{ts}}{\partial z^2} \right|_{\max} = - \left. \frac{\partial F_{ts}}{\partial z} \right|_{\max} > k$$

Only long-range forces measurements (Hamaker constant) are possible




First flexural mode frequency shift



Resonance frequency decreases →
Attractive tip-sample force


Q-factor constant →
Non-dissipative interaction

- Tip-sample not in contact
- Dissipation mechanism probably due to some interaction of the very end of the tip with the surface





$$\frac{\partial F_{ts}}{\partial z} = -2k \frac{\Delta f}{f_0}$$





Non retarded van
der Waals force

$$\frac{\partial F_{ts}}{\partial z} = \frac{HR}{3(z - z_0)^3}$$

$HR = 3 \cdot 10^{-27} \text{ J m}$
 z_0 surface position

H Hamaker constant, R tip radius ($\approx 30 \text{ nm}$), $z - z_0$ tip-sample distance (piezo-tube position, z_0 surface position and cantilever static deflection $\Delta z = F_{ts} / k$)


Second flexural mode



Second mode's stiffness higher than first mode → lower resonance frequency decreases


Flexural mode	1	2	3
k_n (N/m) [1]	0.12	5.0	40
k_n / k_1 (exp.)		43	350
k_n / k_1 (th. [2])		39	300

Non-dissipative interaction




[1] J.E.Sader *et al.*, *Rev. Sci. Ins.* **70**, 3967 (1999)

[2] J.Melcher *et al.*, *Appl. Phys. Lett.* **91**, 053101 (2007)




$$\frac{\partial F_{ts}}{\partial z} = -2k_i \frac{\Delta f_i}{f_{i0}}$$

$i = 1, 2$




Second mode's spring constant k_2

- by Sader method → No matching



$$\frac{\partial F_{ts}}{\partial z} = -2k_i \frac{\Delta f_i}{f_{i0}}$$

$i = 1, 2$



Second mode's spring constant k_2

- by Sader method → No matching
- Stiffer spring constant (tip mass effect) → Good matching

Tip mass loading



1st: no significant effect (1%)

2nd: the node shifts toward the free end
the cantilever effective length reduces
the equivalent stiffness increases

- A **tip mass** that is 10% of the cantilever mass nearly doubles the **second mode** equivalent stiffness. [1]
- Sader method [2] is accurate only for the lower modes.


$k = 0.12 \text{ N/m}$ [2]	k_2 Sader meth.	k'_2 fit to k_{ts}	Th. no tip [1]	Th. tip 10% CL [1]
$k_2 (\text{N/m})$	5.0	9.0		
k_2/k	42	75	40.2	74.9

[1] J.Melcher *et al.*, *Appl. Phys. Lett.* **91**, 053101 (2007)
Y.Sugimoto *et al.*, *Appl. Phys. Lett.* **91**, 093120 (2007)

[2] J.E.Sader *et al.*, *Rev. Sci. Ins.* **70**, 3967 (1999)
MICROSCOPIES A SONDE LOCALE
FORUM 2010

- Introduction
- Free cantilever thermal vibrations
- Probing the tip-sample interaction
 - Frequency shift
 - Potential from Boltzmann distribution
 - Mean-square displacement from power spectral density

Boltzmann probability distribution



- $p(u,z)$ probability of observing the tip at a deflection u from the equilibrium position/tip-surface distance $z \rightarrow$
- number of count at the deflection u divided by the total number of count


Boltzmann probability distribution

$$p(u, z) = p_0 e^{-\frac{V(u, z)}{K_B T}}$$

$V(u,z)$ position dependent total potential

W.F. Heinz *et al.*, *J. Phys. Chem. B* **104**, 622 (2000)

D.O.Koralek *et al.*, *Appl. Phys. Lett.* **76**, 2952 (2000)



Harmonic potential

$$V(u, z) = -K_B T \ln \frac{p(u, z)}{p_0}$$

$$V(u, z) = V_C(u, z) + V_{ts}(u, z)$$

V_c cantilever
potential

V_{ts} tip-sample interaction potential

$$V = \frac{1}{2} k u^2 - \frac{1}{2} \frac{\partial F_{ts}}{\partial z} u^2 = \frac{1}{2} k^* u^2 \quad k^* = k - \frac{\partial F_{ts}}{\partial z}$$

Small oscillation amplitude u ($\partial F_{ts} / \partial z$ constant) \rightarrow harmonic potential

Tip-sample force gradient


$$\frac{\partial F_{ts}}{\partial z} = -\frac{\partial^2 V_{ts}}{\partial s^2}$$

W.F. Heinz *et al.*, *J. Phys. Chem. B* **104**, 622 (2000)

D.O. Koralek *et al.*, *Appl. Phys. Lett.* **76**, 2952 (2000)


$$V_{ts} = V - V_C$$

$V(u,z)$ total potential



Anharmonic contribution
near the surface



$V_{ts}(u,z)$ tip-sample
interaction potential



Interaction potential **parabolic** →
mass-spring approximation correct

Boltzmann probability vs frequency shift

Tip-sample forces gradient



Boltzmann probability distribution

$$\Delta \quad \frac{\partial F_{ts}}{\partial z} = -\frac{\partial^2 V_{ts}}{\partial s^2}$$

Frequency shift method

$$HR = 3 \cdot 10^{-27} \text{ J m}$$

from the frequency shift data fitting

$$\frac{\partial F_{ts}}{\partial z} = \frac{HR}{3(z - z_0)^3}$$

- Introduction
- Free cantilever thermal vibrations
- Probing the tip-sample interaction
 - Frequency shift
 - Potential from Boltzmann distribution
 - Mean-square displacement from power spectral density



Mean-square displacement from PSD

A.P.E. Research
NANOTECHNOLOGY

$$\frac{1}{2} K_B T = \frac{1}{2} k \langle u^2 \rangle$$

Equipartition Theorem

k free cantilever spring constant

$\langle u^2 \rangle$ cantilever mean square deflection



The first mode mean square virtual deflection $\langle u'^2 \rangle_1$ and then the effective spring constant k^* are evaluated at different tip-sample distances z

$$k^* = 0.82 \frac{K_B T}{\langle u'^2 \rangle_1}$$

$$\frac{\partial F_{ts}}{\partial z} = k - k^*$$


H.J. Butt *et al.*, *Nanotechnology* **6**, 1 (1995)

R. Lévy *et al.*, *Nanotechnology* **13**, 33 (2002)


11 March 2010

MICROSCOPIES A SONDE LOCALE
FORUM 2010

Time domain → electronic noise adds to the cantilever virtual deflection



Frequency domain → thermal oscillations isolated from white noise background



$\langle u'_1{}^2 \rangle$ equals the integral of the power spectrum of the thermal fluctuation alone (Lorentzian function) → fit minus white-noise background

Mean-square
displacement from P.S.D.


$$k^* = 0.82 \frac{K_B T}{\langle u_1'^2 \rangle}$$

$$\frac{\partial F_{ts}}{\partial z} = k - k^*$$

Frequency shift method

$$\frac{\partial F_{ts}}{\partial z} = \frac{HR}{3(z - z_0)^3}$$

Tip-sample forces gradient




$HR = 3 \cdot 10^{-27} \text{ J m}$
from the frequency shift data fitting

Thermal oscillations

$$F_{ts} = -\frac{HR}{6(z-z_0)^2}$$


$$HR = 3 \cdot 10^{-27} \text{ J m}$$



Intermittent contact →
Short-range forces

H.Hoelscher *et al.*, Phys. Rev.Lett. **83**, 4780 (1999)

11 March 2010



1. Frequency shift ○
2. Potential △
3. Mean square displacement □

Long-range interaction forces of the order of tens of piconewton

G.Malegori *et al.*, J. Vac. Sci. Techn. in press (2010)






Conclusions

Thermal fluctuations of the AFM cantilever can be exploited to probe the tip-sample interaction

- One temporal trace of the cantilever thermal motion
Three different approaches simultaneously used to measure the tip-sample interaction
- The distance dependence and the magnitude of the force observed are in good agreement with a non retarded van der Waals interaction
- Frequency shift method: mass loading effects on the second eigenmode
- Drawback: due to jump to contact only long-range forces are analyzed

- Investigation of other surfaces
- Measurements in liquids
- Automates the spectroscopy
- Torsional modes analysis



Lateral bending mode (as flexural interchanging cantilever width and thickness) [*]

[*] Espinoza-Beltran *et al.*, *New J. Phys.* **11**, 083034 (2009)



Acknowledgements

A.P.E. Research
NANOTECHNOLOGY

- Vittorio Spreafico, Stefano Prato
APE Research S.r.l., Basovizza (Trieste), Italy
<http://www.aperesearch.com>
- Gabriele Ferrini
Università Cattolica del Sacro Cuore
Dipartimento di Matematica e Fisica,
Via Musei 41, Brescia, Italy
gabriele@dmf.unicatt.it
<http://www.dmf.unicatt.it/elphos>
- Fulvio Parmigiani
Sincrotrone Trieste
Università degli Studi di Trieste
Dipartimento di Fisica, Basovizza, Trieste, Italy



Frequency shift method

The resonance frequency of free cantilever is

$$f_0 = \frac{1}{2\pi} \sqrt{\frac{k}{m^*}}$$

k spring constant, m^* effective mass

The tip-sample interaction F_{ts} changes the resonance frequency.

For small amplitude oscillations

$$k^* = k - \frac{\partial F_{ts}}{\partial z} \quad f = \frac{1}{2\pi} \sqrt{\frac{k^*}{m^*}}$$

k^* effective spring constant

The frequency shift is directly related to the tip-sample force gradient
(exact solution)


$$f = f_0 + \Delta f$$

$$\frac{\partial F_{ts}}{\partial z} = k - k^* = k \left(1 - \frac{f^2}{f_0^2} \right)$$


Y.Martin *et al.*, *J. Appl. Phys.* **61**, 4723 (1987)

F.J.Giessibl, *Rev. Mod. Phys.* **75**, 949 (2003)

Free cantilever's torsional modes



Lateral signal




Material properties

$$\frac{f_n^T}{f^L} = \frac{(2n-1)\pi\sqrt{6}}{\alpha_1^2} \sqrt{\frac{1}{1+\nu} \frac{L}{w}}$$

Geometrical parameters

Boundary conditions

Torsional mode	1	2
f_n (kHz)	250	779
f_n / f_1 (exp.)	23	72
f_n / f_1 (th. [1])	24.71	74.13




Lateral bending mode (as flexural interchanging cantilever width and thickness) [2]

[1] H.J. Van Eysden *et al.*, *J. Appl. Phys.* **101**, 044908 (2007)

[2] Espinoza-Beltran *et al.*, *New J. Phys.* **11**, 083034 (2009)


Tip-sample interaction



Experimental setup


Acquisition system

- Bandwidth: > 1 MHz
- Buffer memory: 128MS
- Sensitivity: 50-200 nm/V



[*] J.L.Hutter, *Langmuir* **21**, 2630 (2005)

Boltzmann probability distribution



- $p(u,z)$ probability of observing the tip at a deflection u from the equilibrium position/tip-surface distance $z \rightarrow$
- number of count at the deflection u divided by the total number of count


Boltzmann probability distribution

$$p(u,z) = p_0 e^{-\frac{V(u,z)}{K_B T}}$$

$V(u,z)$ position dependent total potential

W.F. Heinz *et al.*, *J. Phys. Chem. B* **104**, 622 (2000)

D.O.Koralek *et al.*, *Appl. Phys. Lett.* **76**, 2952 (2000)





Potential from Boltzmann distribution

A.P.E. Research
NANOTECHNOLOGY

Boltzmann probability distribution

$$p(u) = p_0 e^{-\frac{V(u)}{K_B T}}$$

- p_0 normalization constant
- $V(u)$ position dependent potential
- $p(u)$ probability of observing the tip at a deflection u from the equilibrium position/tip-surface distance z

Inverting the relation gives the total potential $V(u)$

$$V(u) = -K_B T \ln \frac{p(u)}{p_0}$$

$$V(u, z) = V_C(u, z) + V_{ts}(u, z)$$

- V_c free cantilever potential
- V_{ts} tip-sample interaction potential

W.F. Heinz *et al.*, *J. Phys. Chem. B* **104**, 622 (2000)

D.O. Koralek *et al.*, *Appl. Phys. Lett.* **76**, 2952 (2000)

$$V(u, z) = V_c(u, z) + V_{ts}(u, z)$$


$$V_c = \frac{1}{2} k u^2$$

- V_c free cantilever potential
- V_{ts} tip-sample interaction potential

Small oscillation amplitude u ($\partial F_{ts} / \partial z$ constant) → harmonic potential

$$V = \frac{1}{2} k^* u^2$$

$$k^* = k - \frac{\partial F_{ts}}{\partial z}$$



Tip-sample force gradient

$$\frac{\partial F_{ts}}{\partial z} = -\frac{\partial^2 V_{ts}}{\partial s^2}$$


W.F. Heinz *et al.*, *J. Phys. Chem. B* **104**, 622 (2000)

D.O. Koralek *et al.*, *Appl. Phys. Lett.* **76**, 2952 (2000)

Cantilever near the surface → F_{ts} interacting force

$$m^* \ddot{u} + \gamma \dot{u} + k u = F_{ts}(z) = F_{ts}(z_0) + \frac{\partial F_{ts}(z_0)}{\partial z} u$$

- The constant force $F_{ts}(z_0)$ only displaces the equilibrium position z_0
- The force derivative $\partial F_{ts} / \partial z$ influences the cantilever oscillations



Small amplitude oscillations → the force gradient $\partial F_{ts} / \partial z$ is constant during the cantilever oscillations u' around the new equilibrium position

$$m^* \ddot{u}' + \gamma \dot{u}' + k^* u' = 0$$


$$k^* = k - \frac{\partial F_{ts}}{\partial z}$$

The tip-sample force gradient is directly related to the effective spring constant k^*

V.L.Mironov *Fundamentals of SPM* (2004)


F.J.Giessibl, *Rev. Mod. Phys.* **75**, 949 (2003)

Tip mass loading



1st: no significant effect

2nd: the node shifts toward the free end,
the cantilever effective length reduces
the equivalent stiffness increases



- A **tip mass** that is 10% of the cantilever mass nearly doubles the **second mode's** equivalent stiffness. [1]
- Sader method [2] is accurate only for the lower modes.


$k = 0.12 \text{ N/m}$ [2]	k_2 Sader meth.	k'_2 fit to k_{ts}	Th. no tip [1]	Th. tip 10% CL [1]
k_2 (N/m)	5.0	9.0		
k_2/k	42	75	40.2	74.9

[1] J.Melcher *et al.*, *Appl. Phys. Lett.* **91**, 053101 (2007)
Y.Sugimoto *et al.*, *Appl. Phys. Lett.* **91**, 093120 (2007)

[2] J.E.Sader *et al.*, *Rev. Sci. Ins.* **70**, 3967 (1999)

Tip mass loading

- The tip mass doesn't influence significantly the equivalent stiffness of the first eigenmode but can have dramatic effect on the equivalent stiffness of the higher eigenmodes.
- A **tip mass** that is 10% of the cantilever mass nearly doubles the **second mode's** equivalent stiffness. [1]
- Sader method [2] is accurate only for the lower modes.



$k = 0.12$ N/m [2]	k_2 Sader meth.	k'_2 fit to k_{ts}	Th. no tip [1]	Th. tip 10% CL [1]
k_2 (N/m)	5.0	9.0		
k_2/k	42	75	40.2	74.9

[1] J.Melcher *et al.*, *Appl. Phys. Lett.* **91**, 053101 (2007)
Y.Sugimoto *et al.*, *Appl. Phys. Lett.* **91**, 093120 (2007)

[2] J.E.Sader *et al.*, *Rev. Sci. Ins.* **70**, 3967 (1999)



A.P.E. Research
NANOTECHNOLOGY